Order Form

Barcode Service for Datamatrix Cryo Rack

(3.) Product order

ordered product.

Please specify the type and quantity of racks to be

barcoded. Please file a separate order form for each

Greiner Bio-One offers the Datamatrix Cryo Rack with a laser written linear barcode on its S1 side and/or a corresponding Datamatrix code on its bottom. Both codes carry the same ID and can be used to identify the rack in sample management and biobanking tasks. Both laser-written features may neither disintegrate nor fall off at temperatures as low as -196 °C.

1. Customer details

Customer No.	Phone
Name	Fax
Institution	E-Mail
Street	Shipping Address
City / Zip Code	
Country	

2. Type of barcode

Please choose whether a linear barcode (Code 128 or Code 39), a Datamatrix code or a combination of both shall be applied. Note that if both coding options are chosen, the sequence of the linear barcode and the Datamatrix code will be identical. Please also specify the sequence start and end as well as the dimensions of the code.

Linear barcode Code 39 with check digit or Code 128 Width of barcode (max. 5 cm) cm Height of barcode (max 1.5 cm) cm In center of S1 side of the Datamatrix Rack Human readable text below barcode (w/o output of check digit)	Side view	 CatNo. 803 277 Rack w/o lid For 1, 2 and 4 ml Cryo.s[™] 20 racks per case cases 	
 Maximum number of characters: 10 Datamatrix Code White Datamatrix code on black background Maximum number of characters: 10 	Bottom side view	 CatNo. 803 202 Rack with low profile lid For 1 and 2 ml Cryo.s™ 20 racks per case cases 	
Sequence start Sequence end Number of characters		Rack with high profile lid For 4 ml Cryo.s™ 15 racks per case cases	
Product samples with barcoding are available on request! Therefore please contact your regional sales representative. (4) Final authorisation Please authorise this order with your signature.			
Place, date F	Released by (signature)		



