


Customer Information

VACUETTE® ID Snap Rings




By using the new VACUETTE® ID Snap Rings, tubes for special analysis purposes can be visually identified at first glance in routine procedures. Throughput times in the laboratory can be eased and accelerated, thus contributing to increased working efficiency.

- Identification of emergency samples
- Identification of special tubes for certain tests
- Identification of special tubes intended for external tests
- Differentiation between similar types of tubes (e.g. K2EDTA or K2EDTA)
- Identification of infectious material

The VACUETTE® ID Snap Rings fit all VACUETTE® 13mm tubes as well as 16mm tubes. Three different neon colours are available: orange, green and pink. There are three hooks on the bottom side of the ring, which provide a firm hold to the tube cap.

Attaching the ring to the cap couldn't be simpler: place the ring with one hook in the middle (ring at an angle) and then push onto cap.



Item n°.	Description	Inner pack	Outer pack
458710	VACUETTE® ID Snap Rings - Starter Kit	6 bags of 500, 2 bags/rack (1000 pcs. of each colour)	3000 pcs.
458711	VACUETTE® ID Snap Rings orange	8 bags of 500 pcs.	4000
458712	VACUETTE® ID Snap Rings green	8 bags of 500 pcs.	4000
458713	VACUETTE® ID Snap Rings pink	8 bags of 500 pcs.	4000