


Ordering Sheet Barcode Service for Cryo.s™

Orderer _____
 Institution _____
 Street _____
 Postal code/City _____

Country _____
 Phone _____
 Fax _____
 E-Mail _____

 For detailed information about "barcode" please ask for our brochure **Barcode-ABC** under (+49) 7022/948-0 or send an E-Mail to info@de.gbo.com.

Please tick off which kind of Cryo.s™ and screw cap colour you need:



Internal thread,
1 ml volume,
self standing



Internal thread,
2 ml volume



Internal thread,
2 ml volume,
self standing



Internal thread,
5 ml volume



External thread,
2 ml volume,
self standing



External thread,
3.5 ml volume,
self standing

Colour screw cap:

- transparent yellow red
 green blue

- Selection of different barcode systems: Code 128; Code 39; 2/5 interleave; EAN; UPC and other
- Recommendation: Code 128C (only numerical code, represents two digits at a time)
Advantage: compact data density and high reliability in reading
- New: Data Matrix (2D barcode) very high data density
Advantage: possibility for clear identification (HIBC) of vial content
- All barcodes with check digit for internal verification and higher reliability in reading
- Printed barcode suitable for preservation in liquid nitrogen (-196 °C), resistant to ethanol (70 %) and DMSO (50 %)
- Minimum order: one original carton (300 pieces for large Cryo.s™ and 500 pieces for small Cryo.s™)

Barcode type: _____

Sequence start: _____

Sequence end: _____

Total number of digits: _____

Total amount of Cryo.s™ with barcodes: _____

Comments: _____

Released by (signature): _____

Date: _____

Greiner Bio-One GmbH
 Maybachstraße 2
 D-72636 Frickenhausen
 Phone: (+49) 70 22/9 48-0
 Fax: (+49) 70 22/9 48-514
 E-Mail: info@de.gbo.com
www.gbo.com/bioscience

